

# TEST TUBES



DEEPALI SCIENTIFIC  
WORKS AND TRADERS

**DEEPALI SCIENTIFIC**  
**SOLUTION FOR A BETTER WORLD**

SHARDA NIWAS, SHIV COLONY, NEAR DR. GAWADE DISPENSARY CHOPDA COURT  
ULHASNAGAR-3, DIST. THANE.

[www.deepaliscientific.com](http://www.deepaliscientific.com)



Test Tubes Are Common Glassware Items Used In Laboratories For Various Purposes, Including Chemical Experiments, Sample Storage, And Qualitative Or Quantitative Analysis. They Are Cylindrical In Shape With A Rounded Bottom And An Open Top, Providing A Versatile And Convenient Vessel For Handling Small Quantities Of Liquids Or Substances.

## OUR ALL TEST TUBES PRODUCTS

TEST TUBE With Rim

3-170ml

TUBE FOR CULTURE COLLECTION

7 X 127mm

TEST TUBE Without Rim

1-170ml

FLAT BOTTOM TUBE Without Rim

12x 75mm

FLAT BOTTOM TUBE With Rim

18 X 138mm

FLAT BOTTOM TUBE With Rim

18 X 55mm

FLAT BOTTOM TUBE With Rim

14 X 110mm